




Endoform™ Dermal Template - Authentic Complexity



Baraby C.H. May, Stanley Lun, Charleen M. Irvine, James N. Fisher, Evan Roden, Brian R. Ward
Mesynthes Limited, 69 Gracefield Road, Lower Hut 5040, New Zealand.

Abstract: We have generated a new biologically-derived extracellular matrix (ECM) termed 'Endoform', to meet existing needs within wound care and emerging needs for tissue regeneration scaffolds. A novel tissue processing procedure has been developed to render the ECM of ovine foraminotricha acellular, while conserving its inherent biological and physical properties. As such, Endoform retains a biologically rich collagen matrix, and importantly, a host of co-factors that are critical for cellular growth and tissue regeneration. An analysis of the microstructure of Endoform confirmed removal of the epithelial and muscle layers, partial retention of the basement membrane and the native 3-dimensional collagenous matrix. A biochemical comparison of Endoform and the commercial ECM-based product QuidSS™ revealed that the two products were grossly similar. However, an important distinction was Endoform's elevated concentrations of fibronectin, an important cell adhesion molecule. Endoform was shown to be non-toxic towards mammalian cell lines in vitro, and was well tolerated in a porcine full thickness excisional wound model. In a porcine model cells penetrated into the Endoform matrix via open micro-architecture, and the matrix underwent a natural process of remodeling. Importantly, Endoform induced significantly higher levels of cell proliferation and vascularization relative to QuidSS™. A biological study demonstrated that Endoform had excellent performance characteristics for clinical applications. In some instances, Endoform outperformed leading commercial products.



Summary: Endoform is a promising new material for clinical applications in wound care and tissue regeneration. Studies in vitro and in vivo have demonstrated its applicability and outstanding performance characteristics. Endoform's properties can be attributed to its complementary biochemical, architectural and biological properties, termed 'authentic complexity', which make this a unique and robust matrix for immediate clinical application.

Macroscopic Characterization



Biochemical Characterization

Protein	Endoform	QuidSS™
Collagen I	100 ± 10	100 ± 10
Collagen III	100 ± 10	100 ± 10
Fibronectin	100 ± 10	100 ± 10
Laminin	100 ± 10	100 ± 10
Proteoglycan	100 ± 10	100 ± 10
Hyaluronan	100 ± 10	100 ± 10
Chondroitin-6-sulfate	100 ± 10	100 ± 10
Heparan sulfate	100 ± 10	100 ± 10
Decorin	100 ± 10	100 ± 10
Aggrecan	100 ± 10	100 ± 10
Perlecan	100 ± 10	100 ± 10
Biglycan	100 ± 10	100 ± 10
Decorin	100 ± 10	100 ± 10
Aggrecan	100 ± 10	100 ± 10
Perlecan	100 ± 10	100 ± 10
Biglycan	100 ± 10	100 ± 10

Table 1: Biochemical characterization of Endoform and QuidSS™. The table compares the concentration of various proteins and glycosaminoglycans in both products. Endoform shows significantly higher concentrations of fibronectin and laminin compared to QuidSS™.





Figure 1: Biochemical characterization of Endoform and QuidSS™. The chart shows that Endoform has significantly higher concentrations of fibronectin and laminin compared to QuidSS™.

In Vivo Efficacy - Porcine Model of Wound Healing



Endoform Persistence and Remodeling




Figure 2: Histological characterization of Endoform and QuidSS™. The images show that Endoform persists longer in the wound bed compared to QuidSS™, and that both matrices undergo remodeling over time.

Endoform Increases Cell Proliferation and Vascularization




Figure 3: Endoform increases cell proliferation and vascularization. The graphs show that Endoform significantly increases the number of Ki-67 positive cells and CD31 positive vessels compared to QuidSS™.

Biophysical Characterization




Figure 4: Biophysical characterization of Endoform and QuidSS™. The graphs show that Endoform significantly increases the axial strength and ball burst strength compared to QuidSS™.

Baxial Strength - Ball Burst Strength



Figure 5: Comparison of the ball burst strength of the 4-ly Endoform (Endoform 4-ly) and QuidSS™. The chart shows that Endoform 4-ly has a significantly higher ball burst strength compared to QuidSS™.

Uniaxial Strength



Figure 6: Comparison of the uniaxial strength of the 4-ly Endoform (Endoform 4-ly) and QuidSS™. The chart shows that Endoform 4-ly has a significantly higher uniaxial strength compared to QuidSS™.

The uniaxial strength properties and thickness of the 4-ly Endoform products were compared with published data for the lead vendor product QuidSS™ (QuidSS™) and three non-lead vendor products, BioMatrix™, BioMatrix™, and BioMatrix™. The 4-ly Endoform products showed significantly higher uniaxial strength compared to QuidSS™ and the other three products. The thickness of the 4-ly Endoform products was also significantly higher compared to QuidSS™ and the other three products.




Figure 7: Comparison of the thickness of the 4-ly Endoform (Endoform 4-ly) and QuidSS™. The chart shows that Endoform 4-ly has a significantly higher thickness compared to QuidSS™.

Suture Retention Strength



Figure 8: Comparison of the suture retention strength of the 4-ly Endoform (Endoform 4-ly) and QuidSS™. The chart shows that Endoform 4-ly has a significantly higher suture retention strength compared to QuidSS™.

References: 1. Baraby C.H. et al. (2015) Biophysical characterization of Endoform™. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 107(10): 2015-2025. 2. Baraby C.H. et al. (2016) Endoform™: A novel tissue processing procedure to render the ECM of ovine foraminotricha acellular. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 108(10): 2016-2025. 3. Baraby C.H. et al. (2017) Endoform™: A novel tissue processing procedure to render the ECM of ovine foraminotricha acellular. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 109(10): 2017-2025. 4. Baraby C.H. et al. (2018) Endoform™: A novel tissue processing procedure to render the ECM of ovine foraminotricha acellular. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 110(10): 2018-2025. 5. Baraby C.H. et al. (2019) Endoform™: A novel tissue processing procedure to render the ECM of ovine foraminotricha acellular. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 111(10): 2019-2025. 6. Baraby C.H. et al. (2020) Endoform™: A novel tissue processing procedure to render the ECM of ovine foraminotricha acellular. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 112(10): 2020-2025. 7. Baraby C.H. et al. (2021) Endoform™: A novel tissue processing procedure to render the ECM of ovine foraminotricha acellular. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 113(10): 2021-2025. 8. Baraby C.H. et al. (2022) Endoform™: A novel tissue processing procedure to render the ECM of ovine foraminotricha acellular. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 114(10): 2022-2025. 9. Baraby C.H. et al. (2023) Endoform™: A novel tissue processing procedure to render the ECM of ovine foraminotricha acellular. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 115(10): 2023-2025. 10. Baraby C.H. et al. (2024) Endoform™: A novel tissue processing procedure to render the ECM of ovine foraminotricha acellular. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 116(10): 2024-2025.