**MedEnvoy**

Prinses Margrietplantsoen 33  
- Suite 123 2595 AM The Hague  
The Netherlands

**MedEnvoy UK Limited**

85, Great Portland Street,  
First Floor, London,  
W1W 7LT  
United Kingdom

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**UK REP**

Easy Medical Device Ltd.  
125 Deansgate  
Manchester, UK M3 2LH

**EC REP**

**S4M EUROPE**  
59, rue Castellion  
01100 Oyonnax - France

**AROA BIOSURGERY LTD.**  
2 Kingsford Smith Place  
Airport Oaks, Auckland 2022  
New Zealand

[www.aroabio.com](http://www.aroabio.com)

**AROA™ Endoform™ Natural** and  
**AROA ECM™** are trademarks of  
Aroa Biosurgery Ltd.

PI.9323.04 | ©March 2023

# Endoform™ Natural

Restorative Bioscaffold

## Instructions For Use

# Endoform™ Natural

Restorative Bioscaffold

## EN Product Description

**Endoform™ Natural** is a 100% w/w intact extracellular matrix (ECM), manufactured from ovine (sheep) forestomach tissue. This advanced wound care device retains the innate biological structure of the native ECM, along with associated macromolecules including elastin, fibronectin, glycosaminoglycans and laminin. When rehydrated with wound exudate or sterile saline, **Endoform Natural** transforms into a soft conforming sheet, which is incorporated into the wound over time.

Not made with natural rubber latex.

Prescription only.

## Intended use

**Endoform Natural** is a sterile, single use ovine forestomach-derived extracellular matrix intended to cover, protect, and provide a moist wound environment.

## Clinical benefit

- Addresses underlying inflammation of chronic healing wounds via broad spectrum wound protease modulation.
- The extracellular matrix provides structure for the patient's own cells to build granulation tissue and supports epithelialization.

## Indications for use

**Endoform Natural** is indicated for the management of wounds including:

- Partial and full thickness wounds
- Pressure ulcers
- Venous ulcers
- Diabetic ulcers
- Chronic vascular ulcers
- Tunneled/undermined wounds
- Surgical wounds (donor sites, grafts, post Moh's surgery, post laser surgery, podiatric, wound dehiscence)
- Traumatic wounds (abrasions, lacerations, second-degree burns, and skin tears)
- Draining wounds

## Contraindications

- **Endoform Natural** is derived from an ovine (sheep) source and should not be used on patients with known sensitivity to ovine (sheep) derived material.
- **Endoform Natural** is not indicated for use on third degree burns.

## Risk statements

### PRECAUTIONS

- Do not apply to wounds with uncontrolled clinical infection, acute inflammation, excessive exudate, or bleeding.
- **Endoform Natural** is terminally sterilized via ethylene oxide. Do not use if the pouch is damaged.
- Discard the device if mishandling has caused possible damage or contamination.
- Single use product. Do not re-sterilize. Discard all unused portions. Device reuse may lead to wound infection.

- Always handle **Endoform Natural** using aseptic technique.
- Do not use past expiry date.

## Storage

**Endoform Natural** should be stored at room temperature in a clean and dry area.

## Directions for use

*These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.*

### 1. Wound Preparation

- a. Prepare the wound bed by cleansing, irrigation and if necessary, debride to remove debris, necrotic tissue or infected tissue.

### 2. Device Application

- a. Select a sheet of the **Endoform Natural** which is slightly larger than the wound and apply aseptically.
- b. **Endoform Natural** can be applied as a whole sheet or trimmed so that it contacts the wound margins. Multiple sheets can be used to cover the entire wound bed.
- c. For ease of handling, apply device by placing the dry material in the wound and rehydrating with exudate or sterile saline. When hydrated **Endoform Natural** transforms into a soft conforming sheet. Ensure that the device conforms to the underlying wound bed.
- d. To protect **Endoform Natural** from adhering to the secondary dressing, consider applying a non-adherent dressing over the device to help protect the tissue while facilitating an optimal moist wound healing environment.
- e. Secure the device using an appropriate secondary dressing. The secondary dressing should be changed according to standard of care, taking into account the level of exudate.
- f. **Endoform Natural** can be used in conjunction with compression therapy, off-loading, and negative pressure wound therapy.

### 3. Reapplication

- a. Duration of treatment and reapplication is determined by the physician and depends upon the wound type and conditions.
- b. The biodegradable **Endoform Natural** is absorbed and incorporated into the wound over time. It is recommended to assess the wound within 72 hours. Reapply when the device has been incorporated into the wound. **Endoform Natural** may last from 1 day and up to 7 days in the wound, depending on the wound type and environment.
- c. Carefully cleanse the wound surface in accordance with established procedures. Do not attempt to remove the off-white to golden gel that forms from residual **Endoform Natural** in the wound as it contains extracellular matrix components that assist in wound healing.
- d. If areas of the device are dry upon inspection, rehydrate with sterile saline and leave in place.
- e. It is not necessary to remove any residual **Endoform Natural** during dressing changes. However, if the product has been overlapped onto the peri-wound area, the remaining loose product that has not incorporated into the wound may be gently removed around the edges if desired.
- f. Change the secondary dressings as needed according to standard of care, taking into account the level of exudate, and when **Endoform Natural** is re-applied.

# Endoform™ Natural


## Restorative Bioscaffold

Users within the EU may access a Summary of Safety and Clinical Performance (SSCP) using Basic UDI-DI 0942190EDTBIOSCAFFOLD5Q through Eudamed (<https://ec.europa.eu/tools/eudamed>) or upon request.

Report any serious incident to Aroa Biosurgery via the local authorized representative (e.g. S4M Europe), as well the applicable Member State competent authority or other regulatory authority, in accordance with local requirements.

### Symbols glossary

Symbols contained in <b>ISO 15223-1 Medical devices</b> – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements are indicated below.		
	5.1.1	Manufacturer
	5.1.2	European Commission Authorized Representative
	-	UK Responsible Person
	-	Swiss Authorised Representative
	-	Conformité Européene
	-	UK Conformity Assessed
	5.1.3	Date of manufacture
	5.1.4	Use by/Expiration date
	5.1.5	Lot number
	5.1.6	Catalog number
	5.1.8	Importer
	5.2.3	Sterilized using ethylene oxide

	5.2.6	Do not resterilize
	5.2.8	Do not use if packaging is damaged
	5.2.11	Single sterile barrier system
	5.3.4	Keep dry
	5.4.2	Do not reuse
	5.4.3	Consult instructions for use
	5.4.8	Contains biological material of animal origin
	5.7.7	Medical device
	-	Number of units
	-	Fenestrated
	-	Non-Fenestrated
	-	High Flow
	-	Disc
Rx Only	-	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician

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