

Proliferative Bioscaffold - US Regulatory Information

Item ID / GTIN	SKU	Layers	Size	US FDA 510(k)	Package Type	Shelf Life (days)	Sterilization Method	Storage Instructions	MRI Compatibility
09421905067096	CM04HA0202	Square	2.5 x 2.5 cm	K200413	Carton	730	Ethylene oxide	Store at room temperature in a dry environment.	MRI Safe
09421905067102	CM04HA0505	Square	5 x 5 cm	K200413	Carton	730	Ethylene oxide	Store at room temperature in a dry environment.	MRI Safe
09421907331461	CM03HA002D	Disc	16 mm	K200413	Carton	730	Ethylene oxide	Store at room temperature in a dry environment.	MRI Safe
09421907331478	CM03HA004S	Square	2 x 2 cm	K200413	Carton	730	Ethylene oxide	Store at room temperature in a dry environment.	MRI Safe
09421907331485	CM03HA006R	Rectangle	2 x 3 cm	K200413	Carton	730	Ethylene oxide	Store at room temperature in a dry environment.	MRI Safe
09421907331522	CM03HA012R	Rectangle	4 x 3 cm	K200413	Carton	730	Ethylene oxide	Store at room temperature in a dry environment.	MRI Safe
09421907331539	CM03HA016S	Square	4 x 4 cm	K200413	Carton	730	Ethylene oxide	Store at room temperature in a dry environment.	MRI Safe
09421907331546	CM03HA025S	Square	5 x 5 cm	K200413	Carton	730	Ethylene oxide	Store at room temperature in a dry environment.	MRI Safe
09421907331577	CM03HA050R	Rectangle	10 x 5 cm	K200413	Carton	730	Ethylene oxide	Store at room temperature in a dry environment.	MRI Safe